sartorius

ePET[®] Electronic Pipette

User Manual

sartorius

ePET[®] Electronic Pipette

English 1

Chinese 23

Contents

1. Intended Use.....	2
2. Product Description.....	2
2.1 ePET® Single-Channel Pipettes.....	2
2.2 ePET® Multichannel Pipettes.....	3
3. Pipette Tips.....	3
4. Unpacking.....	3
4.1 Electrical Specifications.....	3
4.2 Charging the Pipette.....	4
5. Pipette Parts and Materials.....	5
6. Pipette Description.....	6
6.1 Controls.....	6
6.2 Start Button.....	6
6.3 Direction Symbols.....	6
6.4 Display.....	7
6.5 Loading and Ejecting Tips.....	7
6.6 Optional Filters.....	7
7. Programming the Pipette.....	8
7.1 Mode Selection and Mode Recall.....	8
7.2 Setting Speeds.....	8
7.3 Pipetting Mode (P).....	9
7.4 Multi-Dispensing Mode (d).....	10
7.5 Diluting Mode (dd).....	11
7.6 Mixing Mode with Pipetting or Diluting.....	12
8. Pipetting Recommendations.....	13
8.1 Dispensing with Blow-out.....	13
8.2 Dispensing without Blow-out.....	13
8.3 Other Recommendations.....	14
9. Storage.....	14
10. Performance Testing.....	15
10.1 Testing the Performance.....	15
11. Maintenance.....	17
11.1 Cleaning and Decontaminating the Pipette.....	17
11.2 In-house Maintenance.....	17
11.3 Battery Replacement.....	18
12. Troubleshooting.....	19
13. Warranty Information.....	21
14. Recycling Instructions (WEEE).....	21
15. Performance Specifications.....	22

1. Intended Use

The ePET® pipette is intended to be used in liquid handling applications for dispensing liquids in the volume range of 0,2 µl to 5 ml. With ePET® pipettes it is recommended to use Sartorius Optifit Tips or SafetySpace™ Filter Tips to ensure the best compatibility and performance.

This liquid handling device is designed and manufactured to be used as general purpose laboratory equipment. Before operating, it is recommended to read this user manual, which contains useful information, for example, for good pipetting practice.

2. Product Description

Your new Sartorius ePET® electronic pipette offers a cost-effective electronic pipetting alternative by using a direct charging system, which requires no charging stand. The self-calibrating, microprocessor-based system reduces the possibility for human error and instrument contamination, by controlling all piston movements. Also, some models offer replaceable Safe-Cone Filters to help prevent contamination and damage.

Its light weight and ergonomic controls take the effort out of pipetting to help reduce the risk of repetitive strain injuries (RSI) that commonly occur in manual pipetting. All ePET® electronic pipettes operate on the air-displacement principle and use disposable tips.

2.1 ePET® Single-Channel Pipettes

Order No.	Volume Range µl	Increment µl	Tip µl
710521ET	0.2 - 10	0.1	10
710011ET	5 - 100	1	200 , 350
710041ET	50 - 1200	5	1200
710501ET	100 - 5000	50	5000

All pipettes include universal AC-adaptor (EU, UK, US/JPN, AUS, KOR and CHN plugs).

2.2 ePET® Multichannel Pipettes

Order No.	Volume Range μl	Increment μl	Tip μl
710201ET	8-Ch 0.2 - 10	0.1	10
710211ET	8-Ch 5 - 100	1	200, 350
710221ET	8-Ch 25 - 250	5	350
710801ET	8-Ch 50 - 1200	10	1200

3. Pipette Tips

It is recommended that Sartorius Optifit Tips or SafetySpace™ Filter Tips are used with Sartorius pipettes. This ensures the best possible compatibility and thus best possible accuracy and precision, as these tips are designed for the pipette in question.

Optifit Tips and SafetySpace™ Filter Tips are also available as Low Retention Tips.

4. Unpacking

The ePET® pipette package contains the following items:

1. Pipette
2. Grease
3. AC/DC Adapter (acts as a recharging unit)
4. Instructions for use
5. Filters and tweezers (pipettes > 10 μl)
6. ISO8655-6 compliant performance certificate

Make sure that all items are included and that no damage has occurred during shipping.

NOTE! The pipette can only be charged with the original AC/DC Adapter supplied with the pipette.

4.1 Electrical Specifications

Battery

- Rechargeable NiMH battery
- Charging time a maximum of 12 hours for an empty battery

AC/DC Adaptor

- Input voltage and main plug according to local requirements
- Output voltage 9 VDC

4.2. Charging the Pipette

An ON/OFF switch is located at the top of the pipette. This switch protects the battery from discharging while the pipette is not in use and not connected to the charger.

1. Switch the pipette ON (the far left position).
2. Connect the charger to a compatible AC outlet and to the pipette's charger port.

NOTE! The charging light will remain illuminated, when properly connected to the AC/DC charger.

3. If the pipette is new, or the battery charge is low, keep the pipette connected to the charging unit for 12 hours to fully charge the battery before continuing use.
4. Display shows \bar{E} Press the START button twice and the pipette is ready for the default pipetting at maximum volume or for program changes.

NOTE! The pipette will charge in both ON and OFF positions. Leaving the pipette switched ON while charging will retain user settings. Switching it OFF will reset the unit to default settings.

NOTE! If the pipette is left ON for several days, the battery will discharge. As a result the display will turn off, and there will be no response from the keyboard or START button. This is because the charge of the battery is below operating level.


Connect the charger to the pipette's charger port

Switch ON the pipette; the power indicator light comes on.

5. Pipette Parts and Materials

ePET® Single-Channel

Charger port cover (PC)

Finger support (at the back) (ASA/PC)

Keyboard (silicon keys)


Start button (POM)

Handle (ASA/PC)

Power ON/OFF switch

ePET® Multichannel

Tip ejector bar (PETP)

Tip cones
 10-250 µl (sealing: TPU)
 10-250 µl (PVDF)
 50-1200 µl (PPS)

Display (glass)

Keyboard (silicon keys)

Start button (POM)

Handle (ASA/PC)

Tip cone housing (ASA/PC)

6. Pipette Description

The control and the programming of the ePET® electronic pipette are performed using the control keys and display.

6.1. Controls

- Power Switch
- Charger Connection
- Battery Compartment

Control Keys:

- S▶** Changes speed
- M** Moves to next mode
- ▲** Increases volume or speed
- ▼** Decreases volume or speed
- E** Confirms setting changes
- *** Switches Mixing on / off

Start Button

- Triggers all mechanical operations
- Colour-coded for easy identification of pipette volume

Tip Ejector Lever

- Squeeze to remove pipette tip(s)

6.2. Start Button

The START button triggers the aspiration and dispensing operations according to the selected operating mode. Only a quick click is required to operate the button. If the START button is held down, the piston will stop in the lowest position until the button is released. This feature is applicable to all modes of operation, except for multiple dispensing (d).

6.3. Direction Symbols

These symbols indicate the direction in which the piston moves upon pressing the START button. The small RIGHT ARROW in the display means that the next function is to aspirate the liquid. The small LEFT ARROW, in turn, indicates the dispensing function in accordance to the selected operating mode.

6.4. Display

The left side of the display is the status indicator. It informs the user about parameters to be programmed, functions to be performed and the number of dispensings available.

The right side of the display is used for programming and displaying the various volumes needed in different operating modes.

Visible seal when tip is loaded tightly

6.5. Loading and Ejecting Tips

Original Sartorius tips are recommended for use with Sartorius ePET® pipettes. Before fitting a tip make sure that the pipette tip cone is clean. Press the tip onto the cone of the pipette firmly, to ensure an airtight seal. The tip is tightly loaded, when a visible sealing-ring has formed between the tip and the tip cone.

Sartorius ePET® pipettes are designed for simple and light attachment and ejection of the tips. To eject the tip, simply place the pipette over the waste container and squeeze the tip ejector lever.

Squeeze the tip ejector lever, to eject the tip

6.6. Optional Filters

The tip cones of some ePET® electronic pipettes allow the use of a removable Safe-Cone Filter as an option. The filter prevents liquids and vapours from entering the pipette. The filter does not affect the calibration of the pipette. Sartorius filter-tweezers should be used to insert or remove a Safe-Cone Filter to avoid contamination through contact with it.

The use of tweezers to insert or remove a Safe-Cone Filter

Order Information

Order No.	Channels	Volume Range μl	Increment μl	Safe-Cone Filter	
				Standard	Plus
710521ET	1	0.2 - 10	0.1	N/A	N/A
710011ET	1	5 - 100	1	N/A	N/A
710041ET	1	50 - 1200	5	721006	721016
710501ET	1	100 - 5000	50	721006	721016
710201ET	8	0.2 - 10	0.1	N/A	N/A
710211ET	8	5 - 100	1	721014	N/A
710221ET	8	25 - 250	5	721014	N/A
710801ET	8	50 - 1200	10	721006	721016

7. Programming the Pipette

Programming is performed using the six-key keyboard and the LC-display. There are three operating modes offering special functions and varying speeds to choose from.

7.1. Mode Selection and Mode Recall

1. Press **M** repeatedly to view the available modes of the pipette.
2. Press **E** when the desired mode is displayed. The pipette is ready for use in the mode selected.

NOTE: The mode can only be activated when the piston is in its home position (arrow-right is lit), and not during an aspiration or dispensing cycle.

7.2. Setting Speeds

1. Press **S** to display the current aspiration speed.
2. Press **▲** or **▼** until desired aspiration speed is displayed ("5" Fast and "1" Slow).
3. Press **E** to confirm selected speed. Display shows the current dispense speed.
4. Press **▲** or **▼** until desired dispense speed is displayed ("5" being Fast and "1" Slow).
5. Press **E** to confirm speed selection.

NOTE: The speed may not be changed during aspiration or dispensing cycles.

NOTE: The default speed is 3 for all speed settings.

7.3. Pipetting Mode (P)

The pipette performs ordinary pipetting with blow-out.


1. Press **M** to display **P**.
2. Press **E** to confirm the mode change.
3. Select the desired pipetting volume by using **▲** to increase, and **▼** to decrease.

NOTE: While holding down **▲** or **▼**, the volume display starts to scroll faster.

4. Press **E** to confirm selection.
5. Position the tip to aspirate and press the **START** button.
6. Position the tip to dispense. Press the **START** button. The tip is emptied with a blow-out and is ready for the next pipetting.

Step 5 Aspiration
(Pipetting Mode)

Step 6 Dispensing
(Pipetting Mode)

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
In the Multi-Dispensing Mode (d), the pipette performs repetitive dispensings of a selected volume. During this operation, the desired volume, plus the automatically selected excess volume, is aspirated into the tip.

1. Press **M** to display **d**.
2. Press **E** to confirm the mode change.
3. Press **▲** or **▼** until the multi-dispense volume is displayed.
4. Press **E** to confirm selection.
5. Press **▲** or **▼** until the desired number of aliquots is displayed.
6. Press **E** to confirm selection.
7. Position the tip to aspirate and press the START button. The * sign and ARROW LEFT are lit to indicate the reset function.
8. Position the tip to discard priming excess and press the START button.
9. To dispense, position the tip, press the START button and repeat until the cycle is complete.
10. Finally, position the tip to discard any remaining excess and press the START button twice.


NOTE: Upon delivering the programmed volume, it is also possible to continue to aspirate and dispense the same volume without emptying in between. To continue, keep the START button pressed down and within one second the direction of the arrow will change. Keeping the button down, place the tip into the liquid again and the sample is aspirated into the tip by releasing the START button.


Step 7 Aspiration (Multi-Dispensing Mode)

Step 8 Discarding priming excess (Multi-Dispensing Mode)

Step 9 Dispensing repeatedly (Multi-Dispensing Mode)

Step 10 Discarding excess (Multi-Dispensing Mode)


7.5. Diluting Mode (dd)

In the Diluting Mode (dd), two different solutions separated by an air gap are aspirated and then dispensed together with an automatic blow-out. The purpose of the air gap is to prevent contamination while aspirating the second volume, but it will not prevent the two liquids from mixing in the tip.


1. Press **M** to display *dd*.
2. Confirm the mode change by pressing **E**.
3. Select the desired diluent volume (volume 1) by pressing **▲** or **▼**.
4. Confirm by pressing **E**.
5. Press **▲** or **▼** to select the sample volume (volume 2) on the display.
6. Press **E** to confirm selection.
7. Position the tip to aspirate volume 1 and then press the START button.
8. With the tip in the air, press the START button again to aspirate air to form an air gap.
9. Position the tip to aspirate volume 2 and press the START button.
10. Finally, position the tip to dispense, and press the START button.


Step 7 Aspiration (Diluting Mode)

Step 8 Air aspiration for air-gap (Diluting Mode)

Step 9 Aspiration of 2nd fluid (Diluting Mode)

Step 10 Dispensing (Diluting Mode)

7.6. Mixing Mode with Pipetting or Diluting (*)

In the Mixing Mode, the piston is automatically moved up and down to mix the liquid in the delivery vessel. The mixing time is controlled by the START button.

1. Use **M** to select either *P* or *dd* on the display.
2. Press ***** to switch on mixing. The display should show either *P** or *dd**.

NOTE: Pressing ***** will alternately switch mixing on and off.

3. Confirm the mode change by pressing **E**.

For Pipetting:

1. Press **▲** to increase, and **▼** to decrease.
2. Press **E** to confirm selection.
3. Position the tip to aspirate and press START button.
4. Position to dispense and press the START button.

For Diluting:

1. Select the desired diluent volume (volume 1) by using **▲** or **▼**.
2. Confirm selection by pressing **E**.
3. Press **▲** or **▼** to select the sample volume (volume 2).
4. Confirm selection by pressing **E**.
5. Position the tip to aspirate volume 1, and press the START button.
6. With the tip in the air, press the START button to aspirate air for an air gap.
7. Position the tip to aspirate volume 2 and press the START button.
8. Finally, position the tip to dispense and press the START button.

To Mix:

1. Position the tip in the solution, then press and hold the START button. The mixing is done automatically as long as the START button is held down.
2. Position the tip to dispense. Press the START button twice.

NOTE: The mixing is done with about 70% of the total volume.

8. Pipetting Recommendations

By using the various operating modes and special functions, several different liquid handling procedures are possible. Modes P and dd feature an automatic blow-out and others leave the excess liquid in the tip. Follow the recommendations below to ensure optimal performance.


8.1. Dispensing with Blow-out

The P and dd modes have an automatic blow-out function, followed by an immediate return of the piston to the "home" position. To avoid accidental aspiration of the liquid back into the tip, it is recommended that the dispensing is always done above the surface of the liquid.


By holding the START button down, during dispensing, the piston will stop in the lowest position. This allows the tip to be placed against the bottom or the wall of the container. Once the liquid is dispensed, the tip can be removed from the container and the START button can be released.

8.2. Dispensing without Blow-out

The pipette will not perform the blow-out function while using the d mode. Therefore, it is recommended that dispensing is always performed with the tip set against the wall or the bottom of the container. The use of the d mode is especially useful when pipetting small volumes or solutions that have a tendency to foam, or have a high viscosity.

Step 1 Mixing
(Mixing Mode)

Step 2 Dispensing
(Mixing Mode)

8.3. Other Recommendations

- Hold the pipette vertically and place the tip a few millimeters into the liquid when aspirating.
- Pre-rinse the tip by filling and emptying the tip five times. This is important when dispensing liquids which have a viscosity and density different from water or a temperature other than ambient.
- Check that the pipette, tip and liquid are at room temperature.
- Avoid contaminating the tip cone.
- Connect the pipette charger when the pipette is not in use. Switch the pipette OFF when not in use, if it is not connected to the charging unit.
- Never strike the tip cone against a tip tray when mounting tips, as this can damage several internal components.
- Do not drop the pipette or AC/DC-adaptor, as this may cause excessive shock.
- Avoid exposing the unit to extreme temperature changes, humidity and dust.
- Avoid rough handling. Moderate pressure is all that is required when using the keyboard or START-button.
- Avoid leaving the pipette laying on its side with liquid in the tip. The liquid might leak into the pipette's mechanism.
- Always dispense against the inside wall of the receiving vessel. Remove the tip by drawing it up along the inside wall.
- Change the filter on the tip cone regularly (after 50 - 250 pipettings).

9. Storage

When not in use, it is recommended that the pipette is stored on a Sartorius Linear stand, with the ePET® pipette in the ON position and connected to the charger. The green charging light should be lit.

Order No.	Item
725620	Linear Stand for all Sartorius pipettes

Linear Stand

10. Performance Testing

It is recommended, the performance of your ePET® pipettes be checked regularly (e.g. every 3 months), and always after in-house maintenance. However, the user should establish a regular testing routine for their pipettes with regard to accuracy requirements of the application, frequency of use, number of operators using the pipette, nature of the liquid dispensed and the acceptable maximum permissible errors established by the user. (ISO 8655-1.)

10.1. Testing the performance

Performance testing should take place in a draught-free room between 15° and 30°C, kept constant within +/- 0.50°C, and with the relative humidity above 50%. The pipette, tips and the test water should have stood in the test room for a sufficient length of time (at least 2 hours) to reach equilibrium with the room's conditions. Use distilled or de-ionised water (grade 3, ISO 3696). Use an analytical balance with a readability of 0.01 mgs. (ISO 8655-6.)

Weighing

1. Adjust the desired test volume V_S .
2. Carefully fit the tip onto the tip cone.
3. Fill the tip with test water and expel it as waste five times to achieve humidity equilibrium in the dead air volume.
4. Replace the tip. Pre-wet the tip by filling it once with test water and expel it as waste.
5. Aspirate the test water, immersing the tip only 2-3 mm below the surface of the water. Keep the pipette vertical.
6. Withdraw the pipette vertically and touch the tip against the inside wall of the test water container.
7. Pipette the water into the weighing vessel, touching the tip against the inside wall of the vessel just above the liquid surface at an angle of 30° to 45°. Withdraw the pipette by drawing the tip 8-10 mm along the inner wall of the weighing vessel.
8. Read the weight in mgs (m_i).
9. Repeat the test cycle until 10 measurements have been recorded.
10. Convert the recorded masses (m_i) to volumes (V_i)

$$V_i = m_i * Z \quad Z = \text{correction factor}$$

(Table follows)

Z-values ($\mu\text{l}/\text{mg}$):

Temp. ($^{\circ}\text{C}$)	Air Pressure (kPa)			
	95	100	101.3	105
20.0	1.0028	1.0028	1.0029	1.0029
20.5	1.0029	1.0029	1.0030	1.0030
21.0	1.0030	1.0031	1.0031	1.0031
21.5	1.0031	1.0032	1.0032	1.0032
22.0	1.0032	1.0033	1.0033	1.0033
22.5	1.0033	1.0034	1.0034	1.0034
23.0	1.0034	1.0035	1.0035	1.0036
23.5	1.0036	1.0036	1.0036	1.0037

11. Calculate the mean volume (\bar{V}) delivered:

$$\bar{V} = (\sum V_i)/10$$

12. To evaluate conformity, calculate the systematic error e_s of the measurement:

in μl : $e_s = \bar{V} - V_s$ V_s = selected test volume

or in %: $e_s = 100 (\bar{V} - V_s)/V_s$

13. To evaluate conformity, calculate the random error of the measurement:

$$\text{as standard deviation } s = \sqrt{\frac{\sum (V_i - \bar{V})^2}{n - 1}}$$

n = number of measurements (10)

or as a coefficient of variation $CV = 100s/\bar{V}$

14. Compare the systematic error (inaccuracy) and random error (imprecision) with the values in the performance specifications (p. 49) or the specifications of your own laboratory. If the results fall within the specifications, the pipette is ready for use. Otherwise, check both systematic and random errors and, when necessary, proceed to recalibrate (Chapter 9.2).

NOTE: Systematic error (inaccuracy) is the difference between the dispensed volume and the selected test volume. Random error (imprecision) is the scatter of the dispensed volumes around the mean of the dispensed volume. (ISO 8655-1.)

NOTE: Sartorius specifications are achieved in strictly controlled conditions (ISO 8655-6). The user should establish own specifications based on the field of use and the accuracy requirements placed on the pipette (ISO8655-1).

NOTE: This method is based on ISO 8655.

11. Maintenance

To maintain the best results with your ePET® pipette, the cleanliness of every unit should be checked daily. Particular attention should be paid to the tip cone(s).

ePET® pipettes have been designed for easy in-house service. However, Sartorius also provides a complete repair and recalibration service, which includes a service report and performance certificate(s). Please return your pipette to your local representative for repair or performance test/recalibration. Before returning it, please make sure that it is free of any contamination. Please advise our Service Representative of any hazardous materials, which may have been used with your pipette.

NOTE: Check the performance of your ePET® pipette regularly e.g. every 3 months and always after in-house service or maintenance.

11.1. Cleaning and Decontaminating the Pipette

Your ePET® pipette should be checked every day for cleanliness. To clean and decontaminate the outer surface of your pipette use ethanol (70%), isopropanol (60%) or a mild detergent as a cleaning agent.

Gently clean the surface of the pipette with a moistened cloth and then wipe dry. Pay special attention to the tip cone. Change the Safe-Cone Filter if needed (Chapter 6.6.)

11.2. In-house Maintenance

1. Remove the tip ejector collar: Gently twist the tip ejector collar anti-clockwise and slide it off.
2. Using 70% ethanol, 60% isopropanol, or a mild detergent, and a soft lint-free cloth, disinfect the tip ejector collar and the tip cone.
3. Unscrew, anti-clockwise, the tip cone and remove it, exposing the piston. In the event the piston sticks to the tip cone, remove the piston with a pair of tweezers.

4. To avoid scratching the surface of the piston use ethanol (70%), isopropanol (60%) or a mild detergent, and a lint-free tissue when cleaning the piston. Leave the parts to dry.

NOTE: For complete decontamination, place the tip cone, tip ejector collar, piston, O-ring and spring into a beaker containing 70% ethanol, 60% isopropanol, or a mild detergent, and leave for at least 30 minutes. Rinse the parts with distilled water, then dry, preferably with warm air.

5. Apply a thin layer of grease to the piston, using the grease provided. Do not use any other grease. Check that no lint or particles are on the surface of the piston. Avoid excess grease, especially at the bottom of the piston.
6. Reassemble the pipette by screwing the piston and tip cone into their places, replace the tip ejector collar, allowing the ejector handle connection to snap into the notch for the ejector collar. Replace the filter if fitted.
7. Test the pipette by pressing the START-button several times. Test the operation of the tip ejector.

NOTE: Check the performance of your Sartorius ePET® pipette regularly, e.g. every 3 months, and always after in-house service or maintenance.

11.3. Battery Replacement

If the battery is not sufficiently charged for proper operation of the pipette, follow these steps for replacement of the battery.

1. Switch the unit "OFF".
2. Remove the top two screws on the back of the pipette and remove the battery cover.
3. Carefully remove the battery by lifting it straight out of the holder.
4. Install the new NiMH-battery by pressing the positive (+) end against the contact spring at the bottom of the holder.
5. Replace the cover and the screws. Do not overtighten the screws.
6. Dispose of the old battery appropriately.

Battery replacement

12. Troubleshooting

The ePET® electronic pipette has a built-in monitoring program to control the performance of each pipetting action. If the error message ER1 appears on the display, this means the pipette has been unable to perform the attempted action properly. If you receive the error message please do the following:

1. Connect the pipette to the charger for 15 minutes.
2. Clear the error message from the display by pressing .
3. Press START button, which will set the pipette to its home position.

NOTE: As this procedure will empty the tip, it is recommended that you remove the tip before resetting the pipette.

The occasional ER1 situations can be caused by electrical outlets that have been switched off, or if the pipette has been in the OFF position during charging.

Repeated occurrence of the the ER1 message is caused by an internal error; a failure to complete the execution of the pipetting. In which case, it would be necessary to return the pipette to your local Sartorius representative for repair.

Troubles, causes and solutions

Trouble	Possible cause	Solution
Droplets left inside the tip	Unsuitable tip	Use original Sartorius tips
	Non-uniform wetting of the tip plastic	Attach new tip
Leakage or pipetted volume too small	Tip incorrectly attached	Attach firmly
	Unsuitable tip	Use original Sartorius tips
	Foreign particles between tip and tip cone	Clean the tip cone, attach new tip
	Instrument contaminated	Clean and grease O-ring and piston, clean the tip cone (*)
	Insufficient amount of grease on piston and O-ring	Grease accordingly
	O-ring not correctly positioned or damaged	Change the O-ring
	Incorrect operation	Follow instructions carefully
	Calibration altered or unsuitable for the liquid	Recalibrate according to instructions
Push button jammed or moves erratically	Instrument damaged	Send for service
	Piston contaminated	Clean and grease O-ring and piston, clean the tip cone
Pipette blocked, aspirated volume too small	Penetration of solvent vapours	Clean and grease O-ring and piston, clean the tip cone
	Liquid has penetrated tip cone and dried	Clean and grease O-ring and piston, clean the tip cone
Tip ejector jammed or moves erratically	Tip cone and/or ejector collar contaminated	Remove ejector collar, clean with 75% ethanol, isopropanol (60%) or a mild detergent

13. Warranty Information

Sartorius products are meant to be used as described in this user manual.

ePET® electronic pipettes are covered by a warranty for two years against workmanship and defects in materials excluding battery.

ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF THE FAULT IS FOUND TO HAVE BEEN CAUSED BY MISTREATMENT, MISUSE, UNAUTHORIZED MAINTENANCE OR SERVICE OR NEGLIGENCE OF REGULAR MAINTENANCE AND SERVICE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS, OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THIS MANUAL OR WITH TIPS OTHER THAN THE MANUFACTURER'S ORIGINAL ONES.

Each ePET® electronic pipette is tested before shipping by the manufacturer. The Sartorius Quality Assurance Procedure guarantees that the ePET® electronic pipette you have purchased is ready for use. Each ePET® electronic pipette is CE marked, fulfilling the requirements of the EN 55014, 1993/EN 55104, 1995/ ISO 13485:2003 and Directive (98/79 EC).

14. Recycling Instructions (WEEE)

In compliance with European Directive, WEEE (2002/96EC) on waste and reduction of hazardous substances in electrical and electronic equipments, this device must not be recycled as unsorted municipal waste. Instead this device must be collected separately in accordance to the local recycling regulations.

Batteries should be disposed of in accordance with local legal regulations. Do not dispose of batteries with domestic waste.

15. Performance Specifications

The manufacturer's specifications (p. 49) are guaranteed only when the manufacturer's original tips are used. The manufacturer's specifications should be used as guidelines when establishing your own performance specifications in accordance with ISO 8655.

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The status of the information, specifications and illustrations in this manual are subject to change. Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice. All trademarks are Sartorius property unless otherwise stated. Patents granted or pending.
Sartorius Liquid Handling Oy
Sartorius Lab Instruments GmbH & Co. KG