


BIO-RAD



ENGLISH	1 – 14
DEUTSCH	15 – 28
FRANÇAIS	29 – 42
ESPAÑOL	43 – 56

BIO-RAD

BIO-RAD LABORATORIES
PO BOX 210, REGENTS PARK
NSW 2143, AUSTRALIA
ph: +61 (0)29914 2800
fax: +61 (0)29914 2889
www.bio-rad.com

CONTENTS

1 - INTRODUCTION
2 - SETTING THE VOLUME
3 - ASPIRATING AND DISPENSING LIQUIDS
4 - PRE-RINSING
5 - DENSE AND VISCOUS LIQUIDS
6 - BIO-RAD PIPETTE TIPS
7 - RECOMMENDATIONS
8 - RECALIBRATION
9 - TROUBLESHOOTING
10 - PIPETTE KIT AND ACCESSORIES
11 - SPARE PARTS

1 - INTRODUCTION

The **BIO-RAD** pipette is a volumetric instrument designed to measure and transfer liquids precisely and safely. It can measure and transfer volumes from 0.1 μl to 5000 μl depending on the model.

The **BIO-RAD** pipette has a digital counter which displays the volume. The volume is adjusted by turning the black knurled adjustment ring, Fig. 1B, and is continuously adjustable within the volume range for the pipette. The maximum volume of the pipette is shown on the pushbutton and corresponds to the model number, (Fig. 1A).

The volume range from 0.1 μl to 5000 μl is covered by seven pipette models.

Model	Volume range [μl]
BR2	0.1 - 2
BR10	0.5 - 10
BR20	2 - 20
BR100	10 - 100
BR200	20 - 200
BR1000	100 - 1000
BR5000	1000 - 5000

BR2, BR10	Measurement and transfer of micro-volumes, DNA sequencing and enzyme-assay applications.
BR20, BR100, BR200, BR1000	Measurement and transfer of general aqueous solution, acids and bases
BR5000	Measurement and transfer of large volumes

The liquid is aspirated into disposable tips attached to the pipette shaft. Disposable tips ensure maximum safety and eliminate possibility of crosscontamination of the liquid samples.

To protect the user from contaminated tips the **BIO-RAD** pipette is equipped with a built-in tip ejector, (Fig.1D).

The construction of the ejector enables the user to set up the length. The adjustable tip ejector accommodates every variety of tips available on the market. When using narrow tubes, it may be necessary to remove the tip ejector. It is simply removed by pulling down.

Regulating Tip Ejector Length

The regulated "H" spacers, included in the box allow for regulating the length of tip ejector by +1mm or +2mm. An "H0" spacer is inserted on a standard basis. The outside shape of the spacer identifies the size change.

Note that BR5000 does not have a built in tip ejector.

The **BIO-RAD** pipette is a high quality instrument which offers excellent accuracy and precision.

The accuracy and precision (repeatability) of liquid volume depend on the quality of tips used. The values for accuracy and precision given in the table below were obtained using

BIO-RAD tips. Only those tips guarantee correct operation of the pipettes and ensure accuracy and precision of liquid sampling.

Model	Cat no.	Volume [μl]	Accuracy [%]	Precision [%]	Tip
BR2	166-0504-EDU	0.2	± 12	≤ 6.0	10 μl
		1.0	± 2.7	≤ 1.3	
		Max 2.0	± 1.5	≤ 0.7	
BR10	166-0505-EDU	Min 0.5	± 4.0	≤ 4.0	
		5.0	± 1.0	≤ 0.8	
		Max 10.0	± 0.5	≤ 0.4	
BR20	166-0506-EDU	Min 2	± 4.0	≤ 3.0	200 μl
		10	± 1.0	≤ 0.7	
		Max 20	± 0.8	≤ 0.4	
BR100	166-050B-EDU	Min 10	± 1.6	≤ 0.8	
		50	± 1.0	≤ 0.4	
		Max 100	± 0.8	≤ 0.3	
BR200	166-0507-EDU	Min 20	± 1.2	≤ 0.6	
		100	± 1.0	≤ 0.3	
		Max 200	± 0.6	≤ 0.3	
BR1000	166-0508-EDU	200	± 0.8	≤ 0.45	1000 μl
		500	± 0.7	≤ 0.25	
		Max 1000	± 0.6	≤ 0.25	
BR5000	166-0509-EDU	Min 1000	± 0.6	≤ 0.3	5000 μl
		2500	± 0.6	≤ 0.3	
		Max 5000	± 0.5	≤ 0.2	

The accuracy and precision are obtained with **BIO-RAD** tips, using a gravimetric method, performing at least 10 measurements of distilled water at a temperature of $20 \pm 1^\circ\text{C}$, as prescribed by the Deutsches Institut für Normung e.V. To verify the performance of your pipette, request document DIN 12650.

The pipette design enables the user to perform the recalibration process according to the rules presented in section 8.


2 - SETTING THE VOLUME

The counter displays three figures to be read from top to bottom. Additional to the figures on the lower wheel are printed graduations to enable a volume setting in the range increment of each pipette model.

Examples of the meaning of the black and red figures:


Pipettes **BR2**:

Red figures at the bottom = $1/100 \mu\text{l}$
Increment = $0.002 \mu\text{l}$


Pipettes **BR10** and **BR20**:

Red figures at the bottom = $1/10 \mu\text{l}$
Increment = $0.02 \mu\text{l}$

Pipettes **BR100** and **BR200**:

Black figures only = $1 \mu\text{l}$
Increment = $0.2 \mu\text{l}$

Pipettes **BR1000** and **BR5000**:

Red figures at the top - ml

Increment = $2 \mu\text{l}$ $10 \mu\text{l}$

The volume of the pipette is set by turning the black knurled adjustment ring, (Fig. 1 B). To obtain maximum accuracy and reproducibility always change the volume from the higher to the lower value:

- For decreasing the volume turn the adjustment ring slowly to the required setting, making sure not to undergo the mark.
- When increasing the volume turn the adjustment ring until the lower figure wheel comes 1/3 of a turn above the required setting. Now turn the ring slowly backward until the setting reaches the desired value. Be sure not to undergo this mark.

3 - ASPIRATING AND DISPENSING LIQUIDS

Place a tip on the shaft of the pipette. See Section 6 for the appropriate tip. Press the tip on firmly using a slight twisting motion to ensure a positive, airtight seal.

Important: Never aspirate liquids into the **BIO-RAD** pipette without a tip attached.

Aspiration

Press the pushbutton to the first positive stop, (Fig. 2A). Holding the pipette vertically, immerse the tip into the sample liquid. The depth to which the tip is immersed in the sample liquid depends on the model.

Model	Immersion depth (mm)
BR2	≤ 1
BR10	≤ 1
BR20, BR100	2 ÷ 3
BR200, BR1000	2 ÷ 3
BR5000	3 ÷ 6

Release the pushbutton slowly and smoothly to aspirate the sample, (Fig. 2B). Wait one second and then withdraw the tip from the liquid. Wipe any droplets away from the outside of the tip using a medical wipe.

Avoid touching the orifice of the tip.

Dispensing

- Place the end of the tip against the inside wall of the vessel at an angle of 10 to 40 degrees.
- Press the pushbutton smoothly to the first stop, (Fig. 2C). Wait one second.
- Press the pushbutton to the second stop to expel any remaining liquid, (Fig. 2D).
- Keeping the pushbutton depressed to the very end, remove the pipette by drawing the tip against the inside surface of the receiving vessel.
- Release the pushbutton to its starting position, (Fig. 2E).
- Eject the tip by pressing the tip ejector button, (Fig. 2F). Remember to change the tip whenever a different kind of liquid is to be sampled.

4 - PRE-RINSING

When pipetting liquids of higher viscosity or lower surface tension than water (e.g. sera or organic solvents), a film of liquid is formed on the inside wall of the pipette tip. This film can create an error. Since the film remains relatively constant in successive pipetting operations with the same tip, this error can be avoided by forming the film before transferring the first sample. This is done by aspirating a sample and dispensing it back into the same vessel. Since the film is already formed, all of the following samples will have better accuracy and repeatability.

This pre-rinsing operation should be repeated when the volume to be aspirated is changed or when a new tip is used.

5 - DENSE AND VISCOUS LIQUIDS

The **BIO-RAD** pipette specifications of accuracy and precision are based on pipetting distilled water. The handling of liquids with physical qualities of density, viscosity and surface tension differing extremely from water may need a gravimetrically checked compensation of the volume setting. Normally the degree of error resulting from heavy or viscous liquids is negligible if the pipetting is done slowly and carefully. It is most important to give the liquids some time to follow the change of pressure by holding the pipette tip in its position for at least 2 sec. after the aspiration and the blow out stroke.

If in extreme cases this method of operation does not result in accurate values, a compensation could be achieved as follows:

Weigh the liquid pipetted when the **BIO-RAD** pipette is set to the nominal value. Then calculate the set-off from the nominal value:

$$\text{Corr: val.} = 2 \times \text{nom. val.} - \frac{m}{\gamma}$$

m - weight of the sample

γ - density of liquid

Check this operation once again and correct if necessary. Note the corrected value for further pipetting the same kind of liquid.

6 - BIO-RAD PIPETTE TIPS

BIO-RAD tips are made from high performance polypropylene and their quality guarantees the precision and accuracy associated with the **BIO-RAD** pipette. Strict control is maintained throughout the manufacturing process to ensure the highest quality.

The accuracy and precision figures for the **BIO-RAD** pipette are only guaranteed when **BIO-RAD** tips are used. The use of inferior quality tips will seriously degrade the performance of the **BIO-RAD** pipette.

Tips 10

These tips are used for volumes between 0.1 μl and 10 μl . They are used with the BR2 and the BR10 models.

Tips 200

These tips are used for volumes between 2 μl and 200 μl . They are used with the BR20, BR100 and the BR200 models.

Tips 1000

These tips are used for volumes between 100 μl and 1000 μl . They are used with the BR1000 model.

Tips 5000

These tips are used for volumes between 1000 μl and 5000 μl . They are used with the BR5000 model.

7 - RECOMMENDATIONS

Observing the following recommendations will ensure maximal possible accuracy and precision of liquid sampling.

- Make sure to operate the **BIO-RAD** pipette slowly and smoothly.
- The depth of immersion in the sample liquid should be the minimum necessary and should remain constant during aspiration.
- The **BIO-RAD** pipette should be held in a vertical position.
- Change the tip when volume setting is changed or when a different liquid is to be aspirated.

- Change the tip if a droplet remains on the end of the tip from the previous pipetting operation.
- Each new tip should be pre-rinsed with the liquid to be pipetted.
- Liquid should never enter the **BIO-RAD** pipette shaft. To prevent this:
 - Press and release the pushbutton slowly and smoothly.
 - Never turn the pipette upside down.
 - Never lay the pipette on its side when there is liquid in the tip.
- Never force the volume setting beyond its recommended limits.
- When pipetting liquids with temperatures different from the ambient temperature, it is recommended to pre-rinse the tip several times before use.
- Do not pipette liquids with temperatures above 70°C.
- When pipetting acids or corrosive solutions which emit vapours, it is recommended to disassemble the shaft and to rinse the piston and seal with distilled water after finishing the pipetting operation.

8 - RECALIBRATION

BIO-RAD pipettes are calibrated by gravimetric method, using **BIO-RAD** tips and distilled water, at the temperature $20 \pm 1^\circ\text{C}$, according to DIN 12650.

If during pipette operation you find that the accuracy error (the difference between the real aspirated volume and the preset volume) exceeds the permissible value given in the table in section 1, the pipette recalibration procedure should be carried out.

Before starting the recalibration it is necessary to check whether the following requirements have been fulfilled during error determination:

- the ambient temperature, and the temperature of the pipette, tips and water was identical

- the density of the liquid used is close to that of distilled water
- the balance with appropriate sensitivity has been used

Volume checked [μl]	Balance sensitivity [mg]
0.1 - 10	≤ 0.001
10 - 100	≤ 0.01
> 100	≤ 0.1

- mg/ μl conversion factor has been taken into account
- the requirements given in sections 3 and 7 have been fulfilled

If the above conditions are satisfied and the accuracy error for selected volume given in section 1 exceeds the permissible value, the pipette recalibration procedure should be carried out.

The recalibration can be performed within one full turn of the key to the right or to the left only.

Recalibration conditions:

- Ambient temperature and the temperature of the pipette, tips and liquid should be within the range 20-25°C and stabilised during weighing within $\pm 0.5^\circ\text{C}$
- Measurements should be conducted using distilled water
- Balance sensitivity should be suitable for the volume to be controlled

Recalibration procedure:

- Set the dose volume depending on the pipette volume according to the following table:

Model	Range of the pipette volumes [μl]	Preset volume [μl]	Permissible volumes [μl]	Volume change ΔV for full turn of the calibration key [μl] (24 increments)
BR2	0.1 - 2	0.2	0.176 - 0.224	0.06
BR10	0.5 - 10	0.5	0.48 - 0.52	0.33
BR20	2 - 20	2	1.92 - 2.08	0.63
BR100	10 - 100	10	9.84 - 10.16	2.50
BR200	20 - 200	20	19.76 - 20.24	6.30
BR1000	100 - 1000	200	198.4 - 201.6	25.00
BR5000	1000 - 5000	1000	994 - 1006	125.00

- Perform 5 aspirations, weigh each one and calculate the average value of the aspirations
- Calculate average aspirated volume in μl multiplying the average aspiration amount [mg] by the distilled water density coefficient [$\mu\text{l}/\text{mg}$], which depends on temperature and pressure according to the following table:

Temperature [$^\circ\text{C}$]	Pressure [hPa]		
	960	1013	1067
20	1.0028	1.0029	1.0029
21	1.0031	1.0032	1.0032
22	1.0032	1.0033	1.0033
23	1.0035	1.0035	1.0036
24	1.0037	1.0038	1.0038
25	1.0039	1.0040	1.0041

If the average aspirated volume exceeds the permissible value, the following should be done:

- Remove the pipetting pushbutton, (Fig. 4A),
- Holding the volume setting knob to protect it against rotation, insert the calibration key into the cuts of the calibration screw, (Fig. 4B),
- Turn the key clockwise to reduce the aspirated volume, or counter-clockwise to increase the volume. One full turn of the calibration key changes the pipette aspiration volume by the amount given in the table, (Fig. 4C),
- Take out the key and fix the pipetting pushbutton, (Fig. 4D).

Determine the average aspirated volume. The average volume should be within the permissible range given in the table. If the volume exceeds the values stated, the recalibration procedure should be repeated.

In case of pipetting the liquids with physical properties considerably different from those of water, follow the rules given in section 5.

More information on calibration procedures can be obtained from your **BIO-RAD** office.

9 - TROUBLESHOOTING

If the pipette aspirates too little liquid or there is a leak from the tip, check for the following possible causes and perform appropriate correcting measures:

- The shaft nut is loose, (Fig. 3F).

Tighten the shaft nut.

- The sealing surface of the shaft is cracked or scored.

Remove the tip ejector, loosen the shaft nut, and inspect the shaft and the piston assembly. In the models BR2, BR10 and BR20, damage of the shaft may also cause a damaged part of the piston assembly. Replace the damaged parts (see Section 11) and reassemble the pipette, tightening the nut carefully.

- Damage of the piston or seal due to prolonged use with aggressive liquids.

Replace the piston assembly, seal, and O-ring (see Section 11). Rinse the inside of the shaft with distilled water. Lubricate the seal and O-ring with silicon lubricant.

- The pipette is reassembled improperly.

Disassemble the pipette and reassemble it, observing the proper sequence of steps, (Fig. 3).

If there is liquid inside the shaft, clean the pipette as follows:

- Remove the tip ejector, unscrew the shaft nut, then remove the shaft, piston assembly, seal, and O-ring. Rinse the removed parts in distilled water and dry. Reassemble the pipette.

Note: The shaft, shaft nut, and tip ejector may be autoclaved at 121°C for 20 minutes at pressure of 1 bar. The seal and O-ring should not be autoclaved.

The handle and outside of the pipette may be cleaned using a cloth dampened with isopropanol.

If there are air bubbles in the aspirated liquid sample in the tip, proceed as follows:

- Dispense the liquid back to the vessel

- Make sure that the tip is properly and air-tightly attached to the pipette shaft
- Immerse the tip to the proper depth into the liquid
- Aspirate the liquid, releasing the pushbutton more slowly

If the air bubbles still appear in the tip, change the tip.

In case of incorrect pipette aspirations carry out the recalibration procedure according to the guidelines given in section 8.

If the pipette malfunction persists after carrying out the above measures, send the pipette to your **BIO-RAD** service representative.

Before sending the pipette, make sure that it is free from any chemical, radioactive, and microbiological contamination that may pose any danger during transportation and repairing. If possible, please clean the pipette.

10 - PIPETTE KIT AND ACCESSORIES

Pipette kit:

The pipettes are delivered in the kits including:

- Pipette
- Instruction manual
- Calibration key
- Pipette stand
- Ejector regulation spacers
- Tips
- Identification labels

The stand assembly diagram is shown in Fig. 5.

Accessories

Tip Guide (Bulk pack)

Cat. No.	Fits pipette
223-9014	BR2, BR10
223-9035	BR20, BR100, BR200
223-9040	BR1000
6401-B10	BR5000

Tip Guide (Racked)

Cat. No.	Fits pipette
223-9354	BR2, BR10
223-9347	BR20, BR100, BR200
223-9350	BR1000

Stand

Cat. No.	Part description
166-0500	4 place stand

11 - SPARE PARTS

All the spare parts indicated in Fig. 3,4 that is:

A: Pushbutton

B: Calibration knob

C: Shaft

D: Ejector

F: Shaft nut

G: Piston assembly

H: Spacer

I: O-ring

J: Seal

K: Calibration key

can be ordered from your **BIO-RAD** representative (type of pipette and name of the part for this pipette should be specified).

Warning: The replacement of the plunger requires conducting of calibration procedure according to section 8.

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