
Finnpipette® BioControl




Instructions for Use
Bedienungsanleitung
Guide d'utilisation et d'entretien
Instrucciones de uso


Thermo

ELECTRON CORPORATION

Thermo Electron Oy
P.O.Box 100, Fin-01621 Vantaa, Finland
Tel. +358-9-329 100, fax -358-9-3291 0414
www.thermo.com/finnpipette

Display informing modes and parameters


Light indicating recharging when pipette is in the stand

Acceptance key


Mode selection key

Keys to adjust volume and number of steps in stepper mode


5


10


13

11


12

14

12-channel 50µl 2205680
 12-channel 300µl 2205690

50µl 10590360
 300µl 10589300
 1131430

15

Three position pipetting speed selection


Dual stop trigger

Recharge and storage stand

Easy release thumb lever to eject tips

Simple to select pipetting modes and volume adjustment

Autoclavable tip cone module

18

CONTENTS

1	PRODUCT DESCRIPTION
2	GETTING STARTED
3	PIPETTE OPERATION
4	MAINTENANCE AND TROUBLESHOOTING
5	PERFORMANCE AND ORDER INFORMATION

1 PRODUCT DESCRIPTION

The Finnpiquette BioControl is an electronically-assisted pipette with interchangeable tip cone modules for a wide range of liquid handling operations.

The dual-stop trigger that activates the Finnpiquette BioControl uses a natural hand movement, increasing comfort and reducing the risk of repetitive stress injuries.

Thanks to the electric motor and electronic control, pipetting is easy and comfortable, yet still fast and accurate.

Normal and stepper mode pipetting are possible with the same pipette. Forward and reverse pipetting are possible after volume adjustment; no other programming is needed. The dual stop trigger controls the pipetting action.

The pipette works on the air displacement principle (i.e. an air interface is left between the piston and the liquid).

The battery is always charged by rapid recharge technique. If needed, the battery can even be recharged over the lunch hour.

The pipette uses detachable, disposable tips, which are easy to eject with a soft-touch lever action ejector. To insure sterility, the entire tip cone module can be detached and autoclaved.

PACKAGE

The complete Finnpiquette BioControl package contains:

- Finnpiquette with detachable tip cone module(s)
- transformer
- recharge stand
- instructions for use
- calibration software diskette
- tube of grease
- sample Finntips

2 GETTING STARTED

Remove the contents from the package (Fig. 1), and verify that all items listed above are included. Inspect for possible shipping damage. Save the packaging. Make sure that the modules are the desired volume range and that the mains voltage of the transformer is correct.

RECHARGING

The pipette battery may be empty when delivered and must be recharged before initial use. Connect the transformer cord to the outlet on the side of the stand. Plug the wall transformer into an outlet matching the specification on the cover (Fig. 2). **Use only the original BioControl transformer.**

Place the pipette on the stand, noting the light beside the pipette display screen.

A green light indicates that the power is on, and the pipette is ready for use. A red light indicates that rapid recharge has been activated. The initial charge takes about one hour.

To ensure maximum battery capacity, the pipette should be stored on the stand when not in use. The power switches off automatically when the pipette is not in use or on the stand. The pipette display is blank when the power is off.

When you pick up the pipette, the power switches on automatically. If the power is going low, the text **LOW BATT** will appear on screen, indicating that the battery should be recharged. If you continue pipetting, an audible warning will sound.

3 PIPETTE OPERATIONS

The Finnpiquette BioControl has four different operating modes: **RESET**, **PIPETTE**, **STEPPER** and **POWER OFF**.

RESET MODE

The pipette automatically switches to reset mode (on-screen text **CALIBRATE**) whenever the snaplock attaching the tip cone module has been released or the pipette has been recharged from empty. To activate the pipette, close the snaplock, choose the appropriate module (on the display) and press the set key. Press the dual-stop trigger to the second stop. Release the trigger, and the pipette is ready for use (Fig. 6).


Fig.6. On-screen text CALIBRATE

PIPETTE MODE

Pipette is the normal operating mode. Forward and reverse pipetting are possible without programming.

STEPPER MODE

In stepper mode, the volume selected is dispensed repeatedly until the maximum capacity of the tip cone module has been delivered.

POWER OFF MODE

The pipette switches to **POWER OFF** mode after the pipette has been motionless for 10 minutes. A slight movement reactivates the pipette, switching it back on. Volume and mode settings remain as before.

CHANGING THE TIP CONE MODULE

To change the tip cone module, first open the snaplock (about 45 degrees), and then pull out the module. To attach the module to the pipette, insert the module in the pipette and close the snaplock. This activates the **RESET** mode. (To activate the pipette, close the snaplock) Next, choose the appropriate module (channel number and max. volume will be blinking on screen). Accept the previous module by pressing the **SET** key or choose another module by pressing the + or - keys. Confirm your selection with the **SET** key. Then press the dual stop trigger to the second stop, so the piston mechanism inside the handle descends to its lowest position. Release the trigger, and the pipette is ready for use (Fig. 6).

TIP EJECTION

To eliminate the risk of cross-contamination, the Finnpiquette BioControl has a specially designed tip ejector lever. Eject tips by pressing down the lever (Fig. 3). If you are left-handed, just rotate the tip ejector lever to the opposite side of the handle (Fig. 4) and press the lever with your left thumb.

SELECTING PIPETTING MODES

Both **PIPETTE** and **STEPPER** pipetting modes are possible with the Finnpiquette BioControl. **PIPETTE**, the normal aspirate and dispense technique, can be used with or without blow-out. In stepper mode, the pipette repeatedly aspirates and dispenses the selected volume.

To change modes, press the **MODE** key (Fig. 5). The words **PIPETTE** or **STEPPER** will blink on the screen. To confirm your choice, press **SET**. If you do not press **SET** within 20 seconds, the pipette will return to the previous setting. You cannot change modes when the trigger is pressed.

PIPETTE:

Set the volume after selecting and confirming **PIPETTE** mode, (Fig. 7). Press the + or - key to change the volume on the display. If you hold down the + or - key, the figures will scroll on screen. Like the mode text, the volume figures will blink until you confirm your selection with the **SET** key. Once the volume is set, the pipette is ready for use.


Fig.7. Display in PIPETTE mode

STEPPER:

After selecting and confirming **STEPPER** mode, (Fig.8), set the volume using the + and - keys. The unit volume figures will blink until you confirm your choice by pressing **SET**. Once the unit volume is set, the number of repetitions will start to blink on the display. The first option is the maximum number of repetitions. Use the + and - keys to choose the appropriate volume. Press **SET** to confirm your choice.

The maximum number of repetitions depends on the size of the tip cone module and the adjusted unit volume. If the maximum number of repetitions is needed, press **SET** when **MAX STEP** appears on screen (Fig. 9).

PISTON SPEED

Before pipetting, select the piston speed by adjusting the switch on top of the pipette, (Fig. 4). The three pipetting speeds are fast, normal and slow. Fast pipetting is for aspirating and dispensing liquid solutions with a viscosity and surface tension similar to water.

Slower speeds are required for liquids with a high viscosity, low surface tension or a tendency to foam.

PIPETTING TECHNIQUES

Figures 10-12: A = Ready position
B = First stop
C = Second stop

FORWARD AND REVERSE PIPETTING IN PIPETTE MODE

In **PIPETTE** mode, forward and reverse pipetting are possible without further programming. The two stops of the pipetting trigger correspond to the position of the piston, just as in manual pipetting.

FORWARD PIPETTING, (Fig.10)

1. Press the trigger to the first stop.
2. Hold the pipette in a vertical position and dip the tip ends in the liquid. Let the trigger return to the ready position. Withdraw the tips from the liquid, touching them against the edge of the reservoir to remove excess liquid.
3. Place the tips against the wall of the receiving vessel. Dispense the liquid by pressing the pipetting trigger to the second stop. The dispense blow-out action will empty the tips. Remove the tips from the liquid, sliding them along the wall of the vessel.
4. Release the trigger to the ready position.

If necessary, change the pipette tips and continue pipetting.

REVERSE PIPETTING, (Fig.11)

The reverse technique is suitable for dispensing liquids with a high viscosity or tendency to foam. The technique is also recommended for dispensing very small volumes.

1. Press the pipetting trigger to the second stop.
2. Holding the pipette in the vertical position, dip the tip ends in the liquid and release the trigger. This action will fill the tips. Remove the tips from the liquid, touching them against the reservoir to remove excess liquid.
3. Dispense the preset volume by gently depressing the trigger to the first stop. Hold the pipetting trigger at the first stop. The liquid remaining in the tips should not be included in the delivery.
4. The remaining liquid should be discarded by pressing the trigger to the second stop, or aspirated with the next pipetting sequence.


Fig.8. Display in STEPPER mode


Fig.9. Display in STEPPER mode

STEPPER TECHNIQUE, (Fig.12)

In **STEPPER** mode, repeated dispensing of a selected unit volume is possible. There are three phases in stepper pipetting: aspirate, dispense and blow-out.

1. Press the pipetting trigger to the second stop.
2. Holding the pipette in a vertical position, dip the tip ends in the liquid and release the pipetting trigger. This action will fill the tips and start the stepper sequence. Remove the tips from the liquid, touching them against the edge of the reservoir to remove excess liquid.
3. Dispense the preset unit volume by gently pressing the trigger to the first stop. Remove the tips from the liquid, sliding them along the wall of the vessel. Releasing the trigger to the ready position will not move the piston.
4. Activate the next delivery by pressing the trigger to the first stop. The number of strokes left will be displayed on screen.
5. After the last stroke, the words **STEPPER BLOW OUT** will appear on screen, indicating that the dispense sequence has been completed. At this point, only the blow-out volume is left in the tips. Activate blow-out by pressing the trigger to the second stop.
6. Release the trigger, and you are ready to start the next stepper sequence.

CALIBRATION

All Finnpiettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water. Normally, the pipettes do not need adjustment, but they are constructed to permit recalibration and adjustment for liquids of different temperature and viscosity.

DEVICE REQUIREMENTS AND TEST CONDITIONS

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

Volume range readable graduation

under 10µl	0.001mg
10-100µl	0.01mg
above 100µl	0.1mg

Test liquid: Water, distilled or deionized, "grade 3" water conforming ISO 3696. Tests are done in a draft-free room at a constant ($\pm 0.5^{\circ}\text{C}$) temperature of water, pipette and air between 20°C to 25°C . The relative humidity must be above 55%. Especially with volumes under 50 µl the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

CHECKING THE CALIBRATION

The pipette is checked with the maximum volume (nominal volume) and with the lower calibration volume, see Table 1. A new tip is first pre-wetted 3-5 times and a series of ten pipettings are done with both volumes. A pipette is always adjusted for delivery (Ex) of the selected volume. Measuring volumes taken from balance is not allowed. If the calculated results are in the selected limits, the calibration of the pipette is correct. With multichannel models use one of the middle channels to ensure an even load in the module.

Procedure:

1. Do 10 pipettings with the maximum volume.
2. Do 10 pipettings with the lower calibration volume volume.
3. Calculate the accuracy (A) and precision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the results are in the limits of Table 1, then the calibration of the pipette is correct. Otherwise the pipette must be adjusted and checked again.

Single Cannel modules							
Range	Increment	Volume			Precision		
		µl	µl	%	s.d. µl	cv%	
5-40µl	0.1µl	40	± 0.240	± 0.6	0.120	0.3	
		5	± 0.125	± 2.5	0.100	2.0	
40-200µl	1µl	200	± 1.20	± 0.6	0.60	0.3	
		40	± 0.44	± 1.1	0.32	0.8	
200-1000µl	1µl	1000	± 6.0	± 0.6	2.0	0.2	
		200	± 1.8	± 0.9	1.2	0.6	
1-5 ml	0.01ml	5000	± 25.0	± 0.5	10.0	0.2	
		1000	± 10.0	± 1.0	4.0	0.4	
2-10 ml	0.01ml	10000	± 50.0	± 0.5	20.0	0.2	
		2000	± 20.0	± 1.0	6.0	0.3	

Multichannel modules

Range	Channel	Increment	Volume μl	Accuracy μl	Accuracy %	Precision s.d. μl	Precision cv%
0.5-10 μl	8	0.1 μl	10	± 0.24	± 2.4	0.16	1.6
			1	± 0.12	± 12	0.08	8.0
5-50 μl	8,12,16	1 μl	50	± 0.75	± 1.5	0.35	0.7
			5	± 0.25	± 5.0	0.10	2.0
50-300 μl	8,12	1 μl	300	± 3.0	± 1.0	0.90	0.3
			50	± 2.3	± 4.6	0.80	1.5
150-1500 μl	8	1 μl	1500	± 15	± 1.0	4.50	0.3
			300	± 12	± 4.0	7.50	2.5

ADJUSTMENT

Adjustment is done in calibration mode. To change the current calibration factors do following steps:

1. Start the BioControl Calibration Software (included with the pipette).
2. Choose the volume range.
3. Enter the calibration volume results to the fields.
4. Take the pipette and push down

MODE-button and keep it down, then push + and - -buttons as well. When all three buttons are pushed down, release **MODE**-button first, then - and + -buttons. Now **CALIBRATE**-text is blinking, push **SET**-button to accept. Now current HK-factor is blinking followed by PK-factor.

5. Enter the factors to OLD FACTORS-fields.
6. Click "Calculate" button, the new factors are now visible.
7. Take the pipette again and change the factors with + or - buttons, then accept the new value by pushing **SET**-button
8. After adjustment check the calibration according to the instructions.

FORMULAS FOR CALCULATING RESULTS**Conversion of mass to volume**

$$V = (w + e) \times Z$$

V = volume (μl)
 w = weight (mg)
 e = evaporation loss (mg)
 Z = conversion factor for mg/ μl conversion

Evaporation loss can be significant with low volumes. To determine mass loss, dispense water to the weighing vessel, note the reading and start a stopwatch. See how much the reading decreases during 30 seconds (e.g. 6mg = 0.2mg/s). Compare this to the pipetting time from taring to reading. Typically pipetting time might be 10 seconds and the mass loss is 2mg (10s x 0.2mg/s) in this example. If an evaporation trap or lid on the vessel is used the correction of evaporation is usually unnecessary.

The factor Z is for converting the weight of the water to volume at test temperature and pressure. A typical value is 1.0032 $\mu\text{l}/\text{mg}$ at 22°C and 95 kPa. See the conversion table on page 36.

Accuracy (systematic error)

Accuracy is the difference between the dispensed volume and the selected volume of a pipette.

$$A = \bar{V} - V_0$$

A = accuracy
 \bar{V} = mean volume
 V_0 = nominal volume

Accuracy can be expressed as a relative value: $A\% = 100\% \times A / V_0$

Precision (random error)

Precision refers to the repeatability of the pipettings. It is expressed as standard deviation (s) or coefficient of variation (cv)

$$S = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n - 1}}$$

s = standard deviation
 \bar{v} = mean volume
 n = number of measurements
 cv is the relative value of standard deviation.
 $cv = 100\% \times s / \bar{v}$

STERILIZATION

The tip cone module can be sterilized by autoclaving it at 121°C (252 °F) and at 1 bar (15 p.s.i.) for a minimum of 20 minutes. No special preparations are required prior to autoclaving. Steam sterilization bags can be used if needed.

Note: Never autoclave the handle of the FinnpiPETTE BioControl.

Frequent autoclaving of the tip cone module can cause some discoloration. This will not affect the accuracy and precision of the pipette.

4 MAINTENANCE AND TROUBLESHOOTING

When the Finn pipette is not in use, make sure it is safely stored on the recharge stand. The pipette should be checked for dust and dirt on the outer surfaces at the beginning of each workday; pay close attention to the tip cone module in this regard. Use only 70% ethanol as a cleaning solvent.

If the pipette is used daily, it should be checked every three months, following the procedure below:

MULTICHANNEL MODULES 5-50 AND 50-300 µl:

1. Detach the tip cone module by releasing the snaplock and removing the module.
2. Use a screwdriver to remove the four screws in the module cover.
3. Pull out the colour code ring and lift out the module spring.
4. Lift the upper end of the tip ejector bar slightly and pull it out.
5. Open the cover and remove the pistons for cleaning (Fig. 15 12-ch, Fig. 16 8-ch).
6. If needed, replace the seal by lifting the cover ring carefully from its snap joint with a screwdriver.
7. Clean the piston, rings and tip cones with a dry nap-free cloth.
8. Grease the cleaned parts with the lubricant that comes with the pipette.
9. Reassemble the cover ring, spring, support ring and O-ring on the piston; place the assembly into the tip cone; and close the cover ring snap joint.
10. Install the piston bar, pistons and tip cone into the cover A, and close the cover B with the four screws.
11. Install the tip ejector and module spring on the neck of the module. Press the spring under the tip ejector. Close the ejector with the colour code ring.

MULTICHANNEL MODULES 1500 µl:

1. Detach the tip cone module by releasing the snaplock and removing the module.
2. Use a screwdriver to remove the four screws in the module cover.
3. Pull out the colour code ring and lift out the module spring.
4. Lift the upper end of the tip ejector bar slightly and pull it out.
5. Open the cover and remove the pistons and tip cones for cleaning (Fig. 19 8-ch 1500µl).
6. Clean the piston, seals and tip cones with a dry nap-free cloth.
7. Grease the cleaned parts with the lubricant that comes with the pipette.
8. If the piston seal is defect, the whole piston assembly has to be replaced.
9. Reassemble the pistons in to the tip cones.
10. Install the piston bar, pistons and tip cones into the cover A, and close the cover B with the four screws.
11. Install the tip ejector and module spring on the neck of the module. Press the spring under the tip ejector. Close the ejector with the colour code ring.

SINGLE CHANNEL MODULES:

5-40 µl and 40-200 µl modules:

1. Detach the tip cone module by releasing the snaplock and removing the module.
2. Pull out down the tip ejector as far as it goes. Two square holes will be visible on the module body.
3. Use a small flat screwdriver to release the snap joints in the square holes.
4. Pull out the piston.
5. Pull out the other parts (Fig. 17).
6. Clean the piston, rings and tip cones with a dry nap-free cloth.
7. Grease the cleaned parts with the lubricant that comes with the pipette.
8. Assemble the parts in reverse order.
9. Install the locking ring by pressing it with the screwdriver.

200-1000 µl module:

1. Detach the tip cone module by releasing the snaplock and removing the module.
2. Pull out down the tip ejector as far as it goes. Two square holes will be visible on the module body.
3. Use a small flat screwdriver to release the snap joints in the square holes. (Fig. 14).
4. Pull out the piston.
5. The piston assembly cannot be disassembled; only the o-ring can be replaced.
6. Clean the cylinder, piston assembly and tip cones with a dry nap-free cloth.

7. Grease the cleaned parts with lubricant that comes with the pipette.
8. Assemble the parts in reverse order.
9. Install the locking ring by pressing it with the screwdriver.

1-5 ml and 2-10 ml modules:

1. Detach the tip cone module by releasing the latch and removing the module.
2. Pull down the tip ejector as far as it goes. Detach the lower part of the tip ejector. (Snap fit). The upper part can be slid upwards and removed.
3. Press the snap joints of the cylinder to release them. Pull out the cylinder.
4. Pull out the piston.
5. The piston assembly cannot be disassembled; only the o-ring can be replaced
6. Clean the cylinder, piston assembly and tip cone with a dry nap-free cloth.
7. Grease the cleaned parts with the lubricant that comes with the pipette.
8. Assemble the parts in reverse order.

Thermo Electron service department should handle all other service.

CAUTION!

If you send the pipette to us or to our local representative for service, please, ensure that the pipette has been decontaminated before returning.

Also, please note that the postal authorities in your country may restrict the shipment of contaminated material by mail.

TROUBLE SHOOTING

The troubleshooting guide below lists possible problems, causes and solutions.

Problem	Possible cause	Solution
Leakage	Tip incorrectly attached	Attach tips firmly
	Foreign bodies between tip and tip cone	Clean tip cone module; attach new tips
	Foreign bodies between the piston, O-ring and tip cone	Clean and grease O-ring and tip cone. Use grease
	Insufficient grease on tip cone and O-ring	Grease accordingly
Inaccurate dispensing	O-ring damaged	Change the O-ring
	Incorrect operation	Follow instructions carefully
	Tips attached incorrectly	Attach tips firmly
Display is blank	Calibration altered, possibly by misuse	
	Inappropriate calibration. High viscosity liquids may require recalibration.	
CALIBRATE text on the display and the pipette won't operate	Battery is discharged	Mount the pipette in the recharge stand, and make sure that the light beside the display is on
	Power is OFF	Move the pipette slightly
Pipette won't operate	Pipette is in reset mode	Lock the latch and select the module. Press the dual stop trigger to the second stop.
	Tip cone module is improperly attached	Release the latch, attach module firmly to the pipette and lock the latch; press the trigger to the second stop.

TECHNICAL SPECIFICATIONS**Single Channel Modules**

Cat.No	Volume μ l	Increment Δ μ l
2205990	5-40 μ l	0.1 μ l
2206000	40-200 μ l	1 μ l
2206010	200-1000 μ l	1 μ l
2206020	1-5 ml	0.01 ml
2206030	2-10 ml	0.01 ml

Multichannel Modules

Cat.No	Volume μ l	Increment Δ μ l
2206050 (8-ch.)	0.5-10 μ l	0.1 μ l
2205240 (8-ch.)	5-50 μ l	1 μ l
2205250 (8-ch.)	50-300 μ l	1 μ l
2205680 (12-ch.)	5-50 μ l	1 μ l
2205690 (12-ch.)	50-300 μ l	1 μ l
2206060 (8-ch.)	300-1500 μ l	1 μ l
2207030 (16-ch.)	5-50 μ l	1 μ l

ORDERING INFORMATION**Finnpipette BioControl SCP**

Code	Description	Module volume/s
4520X00	FP BioControl 1-ch.	5-40, 40-200, 200-1000 μ l
4520X10	FP BioControl 1-ch.	1-5, 2-10 ml

Finnpipette BioControl MCP

Code	Description	Module volume/s
4531X00	FP BioControl 8-ch.	5-50 μ l
4531X10	FP BioControl 8-ch.	50-300 μ l
4531X20	FP BioControl 12-ch.	5-50 μ l
4531X30	FP BioControl 12-ch.	50-300 μ l
4531X40	FP BioControl 8-ch.	0.5-10 μ l
4531X50	FP BioControl 8-ch.	300-1500 μ l
4531X70	FP BioControl 16-ch.	5-50 μ l

x=0=Eur. 1=US 2=Jap. 3=UK 4=Australia

Battery:

Rechargeable battery. Charging time max. 1 hour for an empty battery.

Transformer:

Input voltage is printed on the transformer. Use only the original transformer supplied with the pipette. Specifications are subject to change without notice.

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Finnpipette® and Finntip® are registered trademarks of Thermo Electron Oy.